



Blood Pressure Cuff Standardisation with Welch Allyn FlexiPort™ Technology

True Standardisation in a Snap.

WelchAllyn®

Advancing Frontline Care™

The one cuff designed for all patient care areas

True cuff standardisation is a snap with Welch Allyn FlexiPort™ blood pressure cuffs

When you standardise your facility on Welch Allyn FlexiPort cuffs, you can use one cuff in virtually any patient care area. Simply equip all blood pressure devices with FlexiPort fittings, assign a cuff to your patients upon admittance to the hospital, and snap the FlexiPort fitting directly into the blood pressure device as the patient moves throughout the hospital, eliminating traditional connector and tube variability.

When you standardise on Welch Allyn FlexiPort Blood Pressure Cuffs:

- Every cuff can be used as one- or two-tube cuff – letting you use one cuff in virtually every room and on virtually every device you have in house
- Part numbers are reduced by up to 60 percent, simplifying ordering and inventory
- Cuffs are available in sizes to fit virtually every patient
- The risk of cross-contamination is reduced by using one cuff per patient
- You get the option of choosing from either disposable or reusable material
- The confusion of multiple tube and connector configurations is eliminated

How can FlexiPort help you?



Clinicians

FlexiPort cuffs make finding the right blood pressure cuff easier and faster than ever before. Spend less time searching for the right cuff, applying it, and connecting it to the device – and have more time to care for your patients.



Biomedical Engineers

Fewer parts and a single connection point for your cuffs mean a more reliable system – and fewer complaints from hospital staff.



Purchasing

Eliminate different tube and connector fittings, reduce part numbers by up to 60%, lower inventory levels, and make the ordering process faster and simpler than ever before.



Infection Control

Reduce the risk of cross-contamination by having one cuff travel with each patient from room to room, floor to floor, and department to department – use one cuff with many devices, not many patients.





With FlexiPort, you can use one cuff on all compatible electronic and manual blood pressure devices – even those that are not manufactured by Welch Allyn.



Why FlexiPort?

Current cuff standardisation methods typically involve a confusing array of cuff types, Y-tubes, clamps, and connectors. This translates into large inventory levels, confusion for the caregiver, and potential device-performance issues due to incorrect or leaky connections. FlexiPort technology provides the only single-point connection solution designed to reduce part numbers and confusion while ensuring reliable, air-tight connections and optimal system performance.





Reduce your part numbers by 60%

Standardising on Welch Allyn FlexiPort blood pressure cuffs makes ordering and inventory of blood pressure cuffs simpler and easier than ever before. Purchase FlexiPort cuffs without tubes or connectors, and your facility will have true cuff standardisation in a snap.



3 simple steps to standardisation

Step 1

Attach FlexiPort fittings to all your blood pressure devices.



Step 2

Choose a cuff material.

Soft Disposable Reusable



Step 3

Select the appropriate cuff size.



Size	Range	Soft Disposable	Reusable
Small Infant (6)	7-10 cm	SOFT-06	REUSE-06
Infant (7)	9-13 cm	SOFT-07	REUSE-07
Small Child (8)	12-16 cm	SOFT-08	REUSE-08
Child (9)	15-21 cm	SOFT-09	REUSE-09
Small Adult (10)	20-26 cm	SOFT-10	REUSE-10
Adult (11)	25-34 cm	SOFT-11	REUSE-11
Adult Long (11L)	25-34 cm	SOFT-11L	REUSE-11L
Large Adult (12)	32-43 cm	SOFT-12	REUSE-12
Large Adult Long (12L)	32-43 cm	SOFT-12L	REUSE-12L
Thigh (13)	40-55 cm	SOFT-13	REUSE-13

Part numbers above represent cuff configurations without tubes and connectors. For a complete listing, contact your local representative.

Neonatal* Soft Disposable Cuffs with Male Luer Slips

Size	Range	One-Tube (10 pack)	Two-Tube (10 pack)
NEO 1	3.3-5.6 cm	5082-101-1	5082-101-2
NEO 2	4.2-7.1 cm	5082-102-1	5082-102-2
NEO 3	5.4-9.1 cm	5082-103-1	5082-103-2
NEO 4	6.9-11.7 cm	5082-104-1	5082-104-2
NEO 5	8.9-15.0 cm	5082-105-1	5082-105-2
Multi-Pack	1 of each	5082-241-9	5082-241-10

*Neonatal cuffs do not have FlexiPort technology.

For more information on how your facility can benefit from FlexiPort cuffs, contact your local Welch Allyn representative or visit us online at www.welchallyn.com/flexiport.

UK Customer Services

Welch Allyn UK Ltd, Cublinton Road, Aston Abbots, Buckinghamshire HP22 4ND, England
Tel: +44 (0) 207 365 6780 Fax: +44 (0) 207 365 9694

Welch Allyn South Africa (Pty) Ltd, 394 Surrey Avenue, Ferndale, Randburg, 2194, South Africa
P O Box 809, Fountainsbleau, South Africa 2032
Tel: +27 (0) 11 777-7555 Fax: +27 (0) 11 777-7556

Welch Allyn Customer Services

Navan Business Park, Dublin Road, Navan, Co. Meath, Republic of Ireland
Tel: +353 46 90 67790 Fax: +353 46 90 67754